

# **NEW VENUS INDUSTRIES**

**111(Old No.47), Kodambakkam Road, Mettupalayam  
Chennai–600 033, Ph: +91–44–23813819/4487/2619**

**E-Mail:** [venusindmng@yahoo.co.in](mailto:venusindmng@yahoo.co.in)



## **PRODUCT PORTFOLIO**

**AUTOCLAVE**



**BALANCE SHIMADZU  
BL SERIES**




| Model   | BL3200HL NEW                                    | BL3200H NEW                      | BL2200H   | BL320H NEW                       | BL220H                            | BL320S              | BL620S              | BL320S               |
|---|---|----------------------------------|---|----------------------------------|-----------------------------------|---------------------|---------------------|----------------------|
| Capacity  | 3200g   | 3200g                            | 2200g   | 320g                             | 220g                              | 3200g               | 620g                | 320g                 |
| Minimum display                                     | 0.01g   | 0.01g                            | 0.01g   | 0.001g                           | 0.001g                            | 0.1g                | 0.01g               | 0.01g                |
| Repeatability(standard deviation)                   | $\sigma \leq 0.01g$                             | $\sigma \leq 0.01g$              | $\sigma \leq 0.01g$                             | $\sigma \leq 0.001g$             | $\sigma \leq 0.001g$              | $\sigma \leq 0.06g$ | $\sigma \leq 0.01g$ | $\sigma \leq 0.006g$ |
| Linearity   | $\pm 0.03g$                                     | $\pm 0.03g$                      | $\pm 0.02g$                                     | $\pm 0.003g$                     | $\pm 0.002g$                      | $\pm 0.1g$          | $\pm 0.01g$         | $\pm 0.01g$          |
| Response time (approx.)                             | 1.0 – 1.2 sec.                                  |                                  |   |                                  | 0.8 – 1.0 sec.                    |                     |                     |                      |
| Ambient temperature                                 | 5 – 40°C / 41 – 104F                            |                                  |   |                                  |                                   |                     |                     |                      |
| Temperature coefficient of sensitivity (10 – 35 °C) | $\pm 5\text{ppm}/^\circ\text{C}$<br>(10 – 30°C) | $\pm 3\text{ppm}/^\circ\text{C}$ | $\pm 5\text{ppm}/^\circ\text{C}$<br>(10 – 30°C) | $\pm 3\text{ppm}/^\circ\text{C}$ | $\pm 10\text{ppm}/^\circ\text{C}$ |                     |                     |                      |
| Pan size (mm)                                       | 160 × 124 (mm)                                  |                                  | 100 × 100 (mm)                                  |                                  | 160 × 124 (mm)                    | 100 × 100 (mm)      |                     |                      |
| Dimensions (approx.)                                | 170W × 240D × 75H mm                            |                                  |   |                                  |                                   |                     |                     |                      |

**ELB SERIES**



**BALL MILL**



**BLOOD MIXER**



**BOD INCUBATOR**



**BOMB CALORIMETER**



**BULK DENSITY APPARATUS**



### Technical Data

|                |                                      |
|----------------|--------------------------------------|
| DISPLAY        | 4 Digit Seven Segment Bright Red LED |
| MAXIMUM COUNTS | 9999                                 |
| KEYBOARD       | 6 Soft Touch Membrane Keys           |
| CYLINDER       | 2 Nos., 50/100 ml Cylinders          |
| DROP HEIGHT    | 25 mm (Approx.)                      |

**CENTRIFUGE REMI C 852**



**R - 8C**

**COLONY COUNTER 361 / 362**

**Technical Data**

| MODEL           | 361                        | 362                         | 363                           |
|-----------------|----------------------------|-----------------------------|-------------------------------|
| DIGITAL DISPLAY | 3 digit, 999 Maximum Count | 4 digit, 9999 Maximum Count | 6 digit, 999999 Maximum Count |
| DISH SIZE       |                            | 110 mm                      |                               |

**1361 / 1362 / 1363**

**Technical Data**

| MODEL      | 1361                        | 1362                        | 1363   |
|------------|-----------------------------|-----------------------------|--|
| DISPLAY    | 3 - Digit Seven Segment LED | 4 - Digit Seven Segment LED | 6 - Digit Seven Segment LED  |
| KEYBOARD   |                             |                             | 6 soft touch membrane keys   |
| CORRECTION |                             |                             | Correction count facility available to cancel the unwanted counts    |
| STORAGE    |                             |                             | 100 Result-count storage retained even if the system is switched off |
| AVERAGE    |                             |                             | Average mode facility available                                      |

**COLOR COMPARATOR**

**COLORIMETER**
**Technical Data**

| MODEL                   | 312/313 DIGITAL        | 321/322 ANALOG   | 1312                   | 1311                   |
|-------------------------|------------------------|--|------------------------|------------------------|
| HIGH STD. GLASS FILTERS | 312 - 8<br>313 - 5     | 321 - 8<br>322 - 5   |                        |                        |
| MINIMUM VOLUME          | 1 ml                   | 1 ml   | 1 ml                   | 1 ml                   |
| DISPLAY                 | 2 ½ digit LED          | Analog Meter   | 16 x 2 line LCD        | 2 ½ digit LED          |
| RANGE                   | 400 to 700 nm          | 400 to 700 nm  | 400 to 700 nm          | 400 to 700 nm          |
| OUTPUT                  | OD : 0 to 1.99         | OD: 0 to 1.99, %T: 0 - 100<br>0-100 % T, 0-1.99 A<br>0-19999 conc. & KFactor |                        | 0-100 % T, 0-1.99 A    |
| RESOLUTION              | OD : 0.01              | OD : 0.01, %T : 1  | OD : 0.01, %T : 1      | OD : 0.01, %T : 1      |
| ACCURACY                | 0.5% FS ±1 Digit       | 2% FS ±1 Digit   | 0.5% FS ±1 Digit       | 0.5% FS ±1 Digit       |
| DECTECTOR               | Photo Cell/Photo Diode | Photo Cell/Photo Diode   | Photo Cell/Photo Diode | Photo Cell/Photo Diode |

**CONDUCTIVITY METER DEEP VISION / TDS 651**
**601 / 611**

**Technical Data**

| MODEL                    | 601 TABLE MODEL  | 611 TABLE MODEL  | 651 TABLE MODEL   | 671 TABLE MODEL                          |
|--------------------------|--|--|---|--|
| DISPLAY                  | 3½ Digit LED   | 3½ Digit LED   | 3½ Digit LED  | 3½ Digit LED                             |
| MEASUREMENT              | Conductivity   | Conductivity   | TDS   | Salinity                                 |
| RANGES                   | 0 To 200 µS/cm<br>0 to 2 mS/cm<br>0 to 20 mS/cm<br>0 to 200 mS/cm<br>0 to 1000 mS/cm | 0 To 200 µS/cm<br>0 to 2 mS/cm<br>0 to 20 mS/cm<br>0 to 200 mS/cm<br>0 to 1000 mS/cm | 0 To 200 ppm<br>0 to 2.000 ppt<br>0 to 20.00 ppt<br>0 to 200.0 ppt<br>0 to 1000 ppt | 0 To 2 ppt<br>0 to 20 ppt<br>0 to 50 ppt |
| ACCURACY                 | ±1% FS±1Digit  | ±1% FS±1Digit  | ±1% FS±1Digit   | ±1% FS±1Digit                            |
| TEMPERATURE COMPENSATION | Auto : 0 to 50° C<br>Manual : 0 to 50° C   | Manual : 0 to 50° C  | Manual : 0 to 50° C   | Manual : 0 to 50° C                      |
| RESOLUTION               | 0.1 µS/cm  | 0.1 µS/cm  | 0.1 ppm   | 0.001 ppt                                |

**EUTECH CON 510**


**Specifications: CyberScan CON 510 Meter**

|                           |   |
|---------------------------|---|
| <b>Conductivity Range</b> | 0 – 20.00, 200.0, 2000 uS/cm;<br>20.00, 200.0 mS/cm                                     |
| Resolution                | 0.01, 0.1, 1 uS/cm; 0.01, 0.1 mS/cm   |
| Accuracy                  | ±1% full scale + 1 digit  |
| <b>TDS Range</b>          | 0 – 10.00, 100.0, 1000 ppm; 10.00, 100.0 ppt<br>max 200 ppt depending on factor setting |
| Resolution                | 0.01, 0.1, 1 ppm; 0.01, 0.1 ppt   |
| Accuracy                  | ±1% full scale + 1 digit  |
| <b>Temperature Range</b>  | -10.0 to +110.0 °C (14.0 to 230 °F)   |



TRUST IN QUALITY

DO METER LUTRON 5509



| <b>SPECIFICATIONS</b> |  |
|-----------------------|--|
| Display               | 13 mm ( 0.5" ) LCD, 3 1/2 digits.                                    |
| Measurement Range     | Dissolved Oxygen - 0 to 20.0 mg/L.                                   |
| Resolution            | 0.1 mg/L..   |
| Accuracy              | $\pm 0.4$ mg/L. (after calibration within $23 \pm 5^\circ\text{C}$ ) |
| Sensor Temp.          | Automatic from 0 to $40^\circ\text{C}$ .                             |
| Compensation          |  |
| Panel adj. knob       | ZERO knob, CAL knob.   |
| Battery               | 006P DC 9V battery ( heavy duty ).                                   |

5510



#### SPECIFICATIONS

|                     |  |
|---------------------|--|
| Circuit             | Custom one-chip of microprocessor LSI circuit.   |
| Display             | Dual function meter's display, 13 mm(0.5"). Super large LCD display with contrast adjustment for best viewing angle. |
| Measurement & Range | Dissolved Oxygen 0 to 20.0 mg/L(liter). Oxygen in Air 0 to 100.0 %.  |
| Temperature         | 0 to $50^\circ\text{C}$ .  |
| Resolution          | Dissolved Oxygen 0.1 mg/L. Oxygen in Air 0.1 % O <sub>2</sub> .  |
| Accuracy (23 5 )    | Temperature 0.1 . Dissolved Oxygen 0.4 mg/L. Oxygen in Air 0.7% O <sub>2</sub> .                                     |


FLAMEPHOTOMETER



#### Technical Data

| MODEL              | 381   | 391  |  |
|--------------------|---|--|--|
| RANGE              | Na : 0-100 ppm<br>K : 0-100 ppm<br>Ca : 15-100 ppm (Optional)<br>Li : 10-100 ppm (Optional) | Serum<br>Na : 0-200 meq/l<br>K : 0-10 meq/l<br>Ca : -<br>Li : 0-2 meq/l (Optional) | Urine<br>0-200 meq/l<br>0-100 meq/l(Optional)<br>0-10 meq/l<br>- |
| READOUT            | 2½ digit bright red LED display   | Dual 2½ digit bright red LED display   |  |
| ACCURACY           | $\pm 2\%$ upto 40 ppm<br>$\pm 5\%$ above 40 ppm   | $\pm 2\%$ upto 100 meq<br>$\pm 5\%$ above 100 meq                                  |  |
| MEASUREMENT SYSTEM | One element at a time   | 2 elements on dual LED display   |  |

GLASS DISTILLATION



#### Technical Data

1382 / 1381

|             | ppm   | General Range<br>meq/l   | Serum Range<br>meq/l  | Urine Range<br>meq/l  |
|-------------|---|--|---|---|
| RANGE       | Na : 0-100<br>K : 0-100<br>Ca : 0-100<br>Li : 0-100                                       | Upto 250 meq/l, 1:100 dil.<br>Upto 250 meq/l, 1:100 dil.<br>Upto 250 meq/l, 1:100 dil.<br>Upto 250 meq/l, 1:100 dil. | 0-200 meq/l, 1:100 dil.<br>0-10 meq/l, 1:100 dil.<br>-<br>0-2 meq/l, 1:5 dil. | 0-250 meq/l, 1:100 dil.<br>0-100 meq/l, 1:100 dil.<br>0-10 meq/l, 1:2 dil.<br>- |
| SENSITIVITY | General Mode : Na:0.5 ppm, K:0.5 ppm, Li:0.5 ppm, Ca:15 ppm In 1382 only                  |  |   |   |
| FILTERS     | Na & K Supplied (Ca & Li are optional at extra cost)                                      |  |   |   |
| RESOLUTION  | 0.1 ppm / meq.  |  |   |   |
| DISPLAY     | 20 x 4 Character alphanumeric LCD display with backlit In 1382 only/ Fro 1382 LED Display |  |   |   |

GERBER CENTRIFUGE



## TRUST IN QUALITY

### GLASSWARES



### HAEMOGLOBINOMETER

1315



### Technical Data

|                  |                   |
|------------------|-------------------|
| MEASURING METHOD | Cyanmethemoglobin |
| MEASURING RANGE  | 0-30 gm/dl of Hb  |
| DILUTION RATIO   | 1 : 251           |
| SAMPLE VOLUME    | 1 ml              |
| DISPLAY          | 3 Digit Display   |
| WAVELENGTH       | 546 nm            |

### HOT AIR OVEN



### SIZE OF INNER CHAMBER

| WxHxD                     | LTS | RATING | SHELVES |
|---------------------------|-----|--------|---------|
| 300x300x300mm 12"x12"x12" | 28  | 750    | 1       |
| 355x355x355mm 14"x14"x14" | 45  | 1000   | 2       |
| 455x455x455mm 18"x18"x18" | 95  | 1500   | 2       |
| 455x605x455mm 18"x24"x18" | 125 | 1750   | 2       |
| 605x605x605mm 24"x24"x24" | 224 | 220    | 3       |
| 605x910x455mm 24"x36"x18  | 252 | 2200   | 3       |
| 605x910x455mm 24"x36"x24  | 336 | 2500   | 3       |

Air circulation fan can be fitted in any of the above sizes at extra cost  
Digital temperature controller

### KARL FISHER TITRATOR



### Technical Data

|  |  |
|--|--|
| MODEL                                    | 1760   |
| MEASURING RANGE                          | 10 ppm to 100% of moisture   |
| MEASURING METHOD                         | Karl Fischer volumetric titration method   |
| K.F. DISPENSING                          | Highly sensitive motorized dispensing system   |
| RESOLUTION OF TITRANT VOLUME             | 0.01 ml  |
| INDICATION OF TITRANT VOLUME             | 0.00 to 99.99 ml on display  |
| MOSITURE CONTENT CALCULATION (AUTOMATIC) | % age : (0.00 to 100%)<br>ppm : (10 to 999999)<br>mg H <sub>2</sub> O : (0.00 to 999.99) |

### WHATMAN FILTER PAPAR



### DURALAB FILTER PAPER



### FILTER HOLDER



### INCUBATOR



| Inner Dimension | Watts |
|-----------------|-------|
| 30 x 30 x 30cm  | 100W  |
| 35 x 35 x 35cm  | 100W  |
| 45 x 45 x 45cm  | 200W  |
| 45 x 45 x 60cm  | 300W  |
| 60 x 60 x 60cm  | 400W  |
| 60 x 60 x 90cm  | 500W  |


### MICROSCOPE



### MOISTURE BALANCE ORDINARY AND DIGITAL



### IMPORTED MODELS



### MUFFLE FURNACE



### NEPHELOMETER / TURBIDITY METER

331 / 341 / 335

| Technical Data |                           |                          |                                 |
|----------------|---------------------------|--------------------------|---------------------------------|
| MODEL          | 331                       | 341                      | 335                             |
| PRODUCT        | Digital Turbidity Meter   | Digital Nephelometer     | Turbidity                       |
| DISPLAY        | 3½ Digital LED            | 3½ Digital LED           | 3½ Digital LED                  |
| DETECTOR       | Photocell/Photodiode      | Photocell/Photodiode     | Photodiode                      |
| RANGE          | 0 to 1000 NTU in 2 ranges | 0 to 200 NTU in 2 ranges | 0-2, 20, 200<br>1000 (4 ranges) |
| RESOLUTION     | 1 NTU                     | 0.01 NTU                 | 0.01 NTU                        |
| ACCURACY       | +3% FS±1digit             | ±2% FS±1digit            | ±2% FS±1digit                   |




## TRUST IN QUALITY

### PETROLEUM PRODUCTS

PENSKY



ABLES



CLEVELAND



CLOUD AND POUR



SMOKEPOINT



PENETRO METER



SOFTINING POINT



### PH METER

111 / 101



#### Technical Data

| MODEL                    | 101   | 111                 |
|--------------------------|---|---------------------|
| DISPLAY                  | 3 ½ Digit LED   | 3 ½ Digit LED       |
| RANGE pH                 | 0 to 14.00  | 0 to 14.00          |
| mV                       | 0 to ±1999  | 0 to ±1999          |
| RESOLUTION pH            | 0.01  | 0.01                |
| mV                       | 1   | 1                   |
| ACCURACY pH              | ± 0.01 pH   | ± 0.01 pH           |
| mV                       | ± 1 mV  | ± 1 mV              |
| TEMPERATURE COMPENSATION | Auto : 0 to 100°C<br>Manual : 0 to 100°C  | Manual : 0 to 100°C |
| ACCESSORIES              | * Combination pH Electrode * Buffer Tablets 4 pH & 7 pH<br>* Buffer Bottles (4 & 7 pH) * Temp. Probe (101E) * Dust Co |                     |

### EUTECH PH METER PH TUTOR



#### METER SPECIFICATIONS

##### CyberScan pH Tutor

|                               |   |
|-------------------------------|---|
| pH Range                      | 0.00 to 14.00 pH  |
| Resolution                    | 0.01 pH   |
| Accuracy                      | ±0.01 pH  |
| Temperature Range             | 0.0 to 100 °C   |
| Resolution                    | 0.1 °C  |
| Accuracy                      | ±0.3 °C   |
| Temperature Compensation      | Automatic / Manual (0 to 100 °C)<br>Default at 25 °C                      |
| Temperature Offset Adjustment | ±5.0 °C   |
| pH Buffer Options             | US Standard (pH 4.01, 7.00, 10.01)<br>NIST Standard (pH 4.01, 6.86, 9.18) |
| No. of Calibration Points     | Up to 3 points  |

pH 510



### POTENIOMETER

118



#### Technical Data

| MODEL         | 118                               |
|---------------|-----------------------------------|
| RANGE         | 0 to ± 199.9 mV<br>0 to ± 1999 mV |
| RESOLUTION    | 0.1 mV & 1 mV                     |
| REPEATABILITY | ±1 mV                             |
| ACCURACY      | ±1 mV ±1 digit                    |

### SOIL PH METER PH 220 S



#### SPECIFICATIONS OF METER


|                   |   |
|-------------------|---|
| Display           | LCD, size : 20 mm x 28 mm.<br><i>* User can reverse the display direction by push button.</i>             |
| Measurement Range | 0 to 14 pH x 0.01 pH  |
| Resolution        | 0.01 pH   |
| Accuracy          | ± 0.07 pH(pH5 to pH 9)<br>± 0.1 pH(pH4 to pH4.9, pH9.1 to pH10)<br>± 0.2 pH(pH1 to pH3.9, pH10.1 to pH13) |

### SEMI AUTO ANALYZER GB 235



### SPECTROPHOTOMETER DOUBLE BEAM

2375



### SOIL MOISTURE METER PMS 714



| SPECIFICATIONS      |   |
|---------------------|---|
| Applications        | Designed to check the moisture level of soil or other similar material.   |
| Measuring Principal | Used the 2 pins electrode to measure the conductive ability of the species, then converter to the reading of " % " Moisture content " of soil sample. |
| Display             | LCD size : 28 mm x 19 mm.   |
| Measurement Range   | 0 % to 50 % moisture content on soil.   |
| Resolution          | 0.1 %.  |
| Accuracy            | ± ( 5 % + 5 d ) F.S.<br>@ 23± 5 °C, F.S. : full scale.  |

6300



#### Technical Data

|                          |  |   |
|--------------------------|--|---|
| MODES                    | 2375 / 6300                            | Basic, Wavelength Scan, Multi Wavelength, DNA/Protein Kinetics, Quantitative, |
| WAVELENGTH RANGE         | 190 - 1100 nm                          |   |
| SPECTRAL BANDWIDTH       | 1.8 nm                                 |   |
| OPTICAL SYSTEM           | Double Beam, Grating 1200 lines/mm     |   |
| WAVELENGTH ACCURACY      | ± 0.3nm                                |   |
| WAVELENGTH REPEATABILITY | ± 0.2 nm                               |   |
| SCANNING SPEED           | HIGH, MED, LOW, MAX 3000 nm/min        |   |
| PHOTOMETRIC RANGE        | 0 - 200 %T, -0.3 to 3A, 0 - 9999 Conc. |   |

SINGLE BEAM 2373



| Technical Data           |  | 2373          | 2371 |
|--------------------------|--|---------------|------|
| WAVELENGTH RANGE         | 190 - 1100 nm                          | 200 - 1100 nm |      |
| SPECTRAL BANDWIDTH       | 2 nm                                   |               |      |
| OPTICAL SYSTEM           | Single Beam, Grating 1200 lines/mm     |               |      |
| WAVELENGTH ACCURACY      | ± 0.5 nm                               |               |      |
| WAVELENGTH REPEATABILITY | 0.3 nm                                 | 0.4 nm        |      |
| WAVELENGTH RESOLUTION    | ± 0.1 nm                               |               |      |
| PHOTOMETRIC ACCURACY     | ± 0.3 %T or ± 0.003A@1A                |               |      |
| PHOTOMETRIC RANGE        | - 0.3 to 3A, 0 - 200 %T, 0 - 9999 Conc |               |      |
| STRAY LIGHT              | 0.05 %T @ 220 nm, 340 nm               | < 0.3 %T      |      |

VISIBLE



2306

| Technical Data            |                                     | 2306          | 2305 |
|---------------------------|-------------------------------------|---------------|------|
| WAVELENGTH RANGE          | 320 - 1000 nm                       | 325 - 1000 nm |      |
| BANDWIDTH                 | 4 nm                                |               |      |
| WAVELENGTH ACCURACY       | ± 1 nm                              | ± 2 nm        |      |
| WAVELENGTH REPEATABILITY  | 0.3 nm                              |               |      |
| WAVELENGTH SETTING        | Auto                                | Manual        |      |
| PHOTOMETRIC ACCURACY      | ± 0.3 %T                            | ± 0.5 %T      |      |
| PHOTOMETRIC REPEATABILITY | 0.2 %T                              | 0.3 %T        |      |
| PHOTOMETRIC RANGE         | - 0.3 - 3.0A, 0 - 200 %T, 0 - 9999C |               |      |

2305



TABLET DISINTEGRATION



| DESCRIPTION    | SDA-273                       | SDA-274                       |
|----------------|-------------------------------|-------------------------------|
| Up-down Speed  | 30+ -2 rpm                    | 30+ -2 rpm                    |
| Up-down Length | 55+ -5mm                      | 55+ -5mm                      |
| Time Display   | 1 Digit For Hr And 2 For Min. | 1 Digit For Hr And 2 For Min. |
| No. Of Baskets | One                           | Two                           |
| Electric Motor | Low rpm A.C. Motor            | Low rpm A.C. Motor            |
| Audio Alarm    | Yes                           | Yes                           |
| Operation Time | Programmable                  | Programmable                  |

TABLET FRIABILITY



| DESCRIPTION            | SFA-263                        | SFA-264                        |
|------------------------|--------------------------------|--------------------------------|
| Rotation Speed         | 25+ -1 RPM                     | 25+ -1 RPM                     |
| Time Display           | 1 Digit For Min And 2 For Sec. | 1 Digit For Min And 2 For Sec. |
| No. Of Drums           | One                            | Two                            |
| Electric Motor         | Low Rpm A.C. Motor             | Low Rpm A.C. Motor             |
| Audio Alarm            | Yes                            | Yes                            |
| Operation Time         | Programmable                   | Programmable                   |
| Auto Removal Of Sample | Yes                            | Yes                            |

TABLET DISSOLUTION

1918



912

ULTRASONIC CLEANER



UV CABINET

| Technical Data  |  | 1912          | 912        |
|-----------------|--|---------------|------------|
| STATIONS        | 8 Station in Model - 1918<br>6 Station in Model - 1916 |               |            |
| DISPLAY         | 20 x 4 Line Alphanumeric LCD                           | 20x4 line LCD | —          |
| SPEED           | 25 - 200 RPM ± 1 RPM                                   | 25 - 160 RPM  | 25-200RPM  |
| RESOLUTION      | 1 RPM  | ±1 RPM        |            |
| RANGE TEMP      | 30 °C - 50 °C  | 20°C to 40°C  | 10 to 50°C |
| REACTION VESSEL | 1000 ml jar with stirrer pedal and basket              | 1000 ml       | 1000 ml    |
| HEATER          | 1 KW   | 1 KV          | 1 KV       |

WATER BATH



WATER AND SOIL TEST KIT

|             | 1160 | 171 | 172 | 161 | 191 | 151 |
|-------------|------|-----|-----|-----|-----|-----|
| pH          | ✓    | ✓   | ✓   | ✓   | ✓   | ✓   |
| Cond.       | ✓    | ✓   | ✓   | ✓   | ✓   | ✓   |
| TDS         | ✓    | ✓   | ✓   | ✓   | ✓   | ✓   |
| DO          | ✓    | ✓   | ✓   | ✓   | ✓   | ✓   |
| mV (ORP)    | ✓    |     |     |     |     |     |
| Salinity    | ✓    | ✓   | ✓   | ✓   | ✓   | ✓   |
| Turbidity   |      |     | ✓   | ✓   | ✓   | ✓   |
| Colorimeter |      | ✓   | ✓   | ✓   | ✓   | ✓   |
| Temp.       | ✓    | ✓   | ✓   | ✓   | ✓   | ✓   |



# Empower Yourself !



From **Lab Instruments** to **Precision Tools**

From **Gadgets** to **Devices**

From **Inventions** to **Solutions**

Find everything from Industry Expert & Renowned Players,

**all at one place!**